

Condenser for lab-based full-field X-ray microscopy with multilayer Laue lenses

S. Niese^{1,2,3}, B. Lechowski¹, D. Rogler¹, R. Dietsch¹, S. Braun⁴, M. Gall³, E. Zschech^{2,3}

1) AXO DRESDEN GmbH, Gasanstaltstr. 8b, 01237 Dresden, Germany; contact@axo-dresden.de

2) Technische Universität Dresden, Dresden Center for Nanoanalysis, Helmholtzstr. 18, 01187 Dresden, Germany

3) Fraunhofer IKTS-MD, Maria-Reiche-Str. 2, 01109 Dresden, Germany

4) Fraunhofer IWS Dresden, Winterbergstr. 28, 01277 Dresden, Germany

Multilayer Laue lenses – MLLs

- Concept: linear zone plate to realize enhanced efficiency and resolution for hard X-ray imaging [1]
- Crossed MLLs for point focusing and full-field imaging
- Fabrication of partial MLLs
 - magnetron sputtering: thin film deposition
 - shaping: precision sawing and focused ion beam milling
 - glue bonding: crossed MLL
 - focal length adjustment by exploiting radial thickness variations
- Properties
 - focal length: 8.0 mm (Cu-K α)
 - zones: 50...2500 = 79...11 nm
 - WSi₂/Si multilayer, height: 48 μ m




Fig 1: Illustration and SEM image of crossed partial multilayer Laue lenses

Full-field X-ray microscopy

- Full-field imaging demonstrated in a laboratory X-ray microscope Zeiss Xradia NanoXCT-100 using the hollow cone illumination of a capillary condenser and one-sided imaging [2]
- Current issues
 - ghost image formation due to Cu-K β radiation
 - reduced intensity due to hollow cone illumination
 - size of the illumination exceeds the field of view





Fig 2: Recent results compared to imaging with a Fresnel zone plate (FZP). A Nickel filter is mandatory to suppress Cu-K β ghost images. [2]

References


- [1] J. Maser, G. B. Stephenson, S. Vogt, W. Yun, A. Macrander, H. C. Kang, C. Liu, R. Conley: Multilayer Laue lenses as high-resolution x-ray optics. Proc. SPIE vol. 5539, p. 185, 2004.
- [2] S. Niese, P. Krüger, A. Kubec, S. Braun, J. Patommel, C. G. Schroer, A. Leson, E. Zschech: Full-field X-ray microscopy with crossed partial multilayer Laue Lenses. Optics Express, 2014. Accepted.

Illumination with AXO ASTIX-f optics

- Two-dimensional focusing X-ray mirror in a modified Montel geometry
 - solid illumination
 - efficient suppression of Cu-K β radiation and bremsstrahlung
 - matching of numerical apertures using a slit system
 - spot simulation considering slope errors
- Properties of the fabricated mirror
 - $f_1 = 325$ mm, $f_2 = 125$ mm
 - designed for Cu-K α radiation
 - secondary focus FWHM = 28 μ m



a) simulation



b) measurement




Fig. 4: Optical path illustration

Experimental set-up

- Demonstrator setup with combined optics for basic evaluation
 - Cu microfocus X-ray tube
 - ASTIX-f mirror & slit system
 - 6-axis sample stage
 - Multilayer Laue lens
 - X-ray detector, 6.45 μ m pixel size
- In progress
 - proof-of-principle demonstration for X-ray microscopy and tomography
 - study of the influence of the illumination
- Outlook — hard X-ray microscopy with ...
 - enhanced efficiency
 - from sub-100 nm to sub-10 nm resolution
 - photon energies $E > 10$ keV



Fig 4: Demonstrator setup – μ Focus Cu X-ray source, ASTIX-f mirror, object, MLL, X-ray detector




Fig. 5: Preliminary results

Acknowledgment

This work is supported by the German Ministry for Education and Research, BMBF, under the Program "IKT 2020 - Research for Innovations", Project No. 16ES0070. This work is partly supported by the German Research Foundation (DFG) within the Cluster of Excellence "Center for Advancing Electronics Dresden" and the German Federal Ministry of Research and Education (BMBF) within the MaKiZu project. Motorized stages are kindly provided by HUBER Diffraktionstechnik GmbH & Co. KG.