

IWS



# Multilayer Laue Lenses with Long Working Distances

Sven Niese<sup>1</sup>, Saeed Mirzaei<sup>2</sup>, Peter Gawlitza<sup>2</sup>

<sup>1</sup> AXO DRESDEN GmbH, Gasanstaltstraße 8b, 01237 Dresden, Germany; sven.niese@axo-dresden.de <sup>2</sup> Fraunhofer IWS Dresden, Winterbergstraße 28, 01277 Dresden, Germany

# Introduction

Multilayer Laue Lenses (MLLs) are innovative diffractive X-ray optics with a high numerical aperture (NA) that enable X-ray focusing to sub-10 nm spatial resolution even for hard x-rays [1]. The principle of MLLs relies on diffraction. The individual layers (zones) are made in such a way, that constructive interference occurs at the focus.

#### Synthesis of Multilayer Laue Lenses

Alternating thin layers must be deposited with a thickness gradient so that they follow zone plate law. Free-standing coatings are structured by short-pulse laser followed by FIB thinning for the target photon energy [2]. Fig. 1 shows the schematic of an MLL.



Fig. 1: Design, fabrication and focusing of Multilayer Laue Lenses.

## Results

#### Application of MLL in synchrotron beamlines

Novel synchrotron X-ray sources provide a very intense beam with extraordinary high brilliance and coherence. The performance of MLLs in focusing a synchrotron beam to a small spot size has been well documented [3,4]. Currently, one major challenge is a limited working distance (WD) impeding the investigation of real specimens. Therefore, the development of MLLs with large WD and acceptable resolution is crucial. Table 1 lists three types of our MLLs with long, very long and extra-long WD (LWD, VLWD, ELWD) developed at Fraunhofer IWS in cooperation with AXO DRESDEN.

Table 2: MLL parameters and properties.

MLL design	Tested		Under test
	LWD	VLWD	ELWD
Materials	Mo/C/Si/C	Mo/C/Si/C	Mo/C/Si/C
Focal length (mm) @ 12 keV	9	45	67.5
Working distance (mm) @ 12 keV	3.1	25	39
Geometry	Flat tilted	Flat tilted	Flat tilted
Stack height (µm)	50	100	150
Individual layers	12000	17000	24000
Focus size FWHM (nm)	22×22	80×80	85×85

An integrated MLL mount was developed: This mount integrates the vertically and the horizontally focusing MLL to a 2D operating lens device and provides all required degrees of freedom for the internal alignment. This significantly reduces the complexity of the final integration to the beamline. Fig. 2 shows all MLL designs installed on the ESRF ID13 beamline and corresponding sample at the focal position. An order sorting aperture blocks the higher-order diffractions and the direct beam.



Fig. 2: Comparison of MLL designs at 12 keV.

## References

- [1] Bajt, Saša, et al., Light: Science & Applications 7.3 (2018): 17162-17162.
- [2] Kubec, Adam et al., Applied Physics Letters 110 (2017): 111905
- [3] Yan, Hanfei, et al., X-Ray Optics and Instrumentation 2010.1 (2010): 401854.
- [4] Kubec, Adam, et al., Journal of Synchrotron Radiation 24.2 (2017): 413-421.
- [5] Susanne Hönig et al., Optics Express 19 (2011): 16324-16329.
- [6] Björn Enders et al., Proceedings of the Royal Society A vol. 472 (2016): 20160640.
- [7] Meindlhumer, Michael, et al., Journal of Vacuum Science & Techn. A 42.2 (2024).

The focusing characteristics is typically evaluated with ptychography [5,6]. Our MLLs were tested at PETRA III, P03, operating at 12.7 keV. These lenses exhibited the focus size of 30×40 nm<sup>2</sup> and as shown in Fig. 4, the minimal features of reconstructed test patterns could be resolved. Fig. 5 also shows preliminary results from experiments at ID13, ESRF.



Fig. 4: Experimental setup and ptychographic reconstruction at PETRA III, P03.



Fig. 5: MLLs developed at IWS exhibited ≈ 22 nm (LWD design, left) and 53 nm (VLWD design, right) focus size at ESRF, ID13 (preliminary result).

#### MLLs for structural analysis

VLWD MLLs were used to focus the synchrotron beam at ID13, ESRF and investigate stress and structural properties in transistors. A cross-sectional nano-XRD setup as insitu micromechanical experiments is shown in Fig. 6 (a), where a sample is scanned with an MLL focused X-ray beam with ≈ 75 nm FWHM [7]. A notched cantilever with Pt dots as markers is shown in Fig. 6 (b). Local defect density and stress distributions were investigated under load, see in Fig. 6 (c)



Fig. 6: Results from in-situ micromechanical experiments [7].

## Conclusion

- MLLs address the limits of traditional X-ray lenses such as limited resolutions and diffraction efficiencies particularly for hard x-rays
- MLLs enable high-resolution X-ray studies with nm-scale focusing
- Three designs with increasing working distance were successfully developed
- MLLs can be tailored to specific applications offering flexibility in terms of energy range, focal length and working distances
- The integrated MLL mount allows an easy and safe integration to the beamline

### Acknowledgements

Authors acknowledge Dr. Manfred Burghammer from ESRF, Grenoble, France for his assistance with measurements and insightful contribution to discussion. This research was supported by MI-1216, MI-1855 and MI-1463 long-term projects at ESRF Grenoble, France as well as PETRA III LTP: II-20180007 EC, DESY, Hamburg, Germany.